


The medical utility of laboratory tests is limited by variation arising from any source, including biological variation as well as analytical variation (measurement error).



μ = uncertainty, U = expanded uncertainty = $k * \mu$ using a $k = 2$ for approximately 95% confidence and $k=3$ for approximately 99% confidence
 $X \pm 2 \mu$ includes a true value of X with approximately 95% confidence; between the values $(X-2\mu, X+2\mu) = (X-U, X+U)$

At a level near X , MU is $\pm U$; (95% CI = $X-U$ to $X+U$) The number of significant digits given for a MU should be the same as that used for reported results.
 $X \pm 3\mu$ includes the true value with approximately 99% level of confidence.

Reference change values (RCVs) - Determines whether the difference between two results is negligible due to uncertainty or significant due to a genuine change in the condition of the patient.

$RCV_{in \%} > 2.77 * (\%CV_{Analytical})$ at a 95% Confidence Level

$RCV_{in units} > 2.77 * (sd_{Analytical}) = 2.77 ([\%CV_{Analytical} * test\ result]/100\%) = 2.77 ([\mu * test\ result]/100\%)$

If biological variation (CV_1) is known, then $RCV_{in \%} > 2.77 * \sqrt{(\%CV_{Analytical}^2 + \%CV_1^2)}$ at a 95% Confidence Level

$\%CV_{Analytical} = \mu$ $sd_{Analytical} = \mu_{in\ units}$ $sd_{Analytical} = (\%CV_{Analytical} * mean) / 100\%$

- If the method has a quantitation step, such as an absorbance value for determining a cutoff, measurement of uncertainty must be calculated.
- When calculating combined uncertainties for parameters that are calculated using addition and subtraction e.g. Anion gap, the SD or μ value can be used.
- Similarly when calculating combined uncertainties for parameters that are calculated using division and multiplication e.g. creatinine clearance, the sd or μ must first be converted to %CV.